

## Exact Volume Transfer Pipets

*Easy to use, economical and accurate!*

These **exact volume transfer pipets** are produced from chemically resistant, low density polyethylene and are designed to dispense precise volumes. They are available in the 5 popular sizes shown below. Exact volume pipets are packaged 500 pipets per bag.

### Exact Volume

Actual Size

#### 50 $\mu\text{L}$ (0.05 mL) capacity



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
139110	Non-sterile	500 pcs/bag, 10 bags/case	5000	<ul style="list-style-type: none"><li>Capacity: 50 <math>\mu\text{L}</math></li><li>Bulb draw: 50 <math>\mu\text{L}</math></li><li>Total length: 59 mm (2.3")</li><li>Stem diameter: 2.2 mm (0.1")</li></ul>

#### 75 $\mu\text{L}$ (0.075 mL) capacity



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
139112	Non-sterile	500 pcs/bag, 10 bags/case	5000	<ul style="list-style-type: none"><li>Capacity: 75 <math>\mu\text{L}</math></li><li>Bulb draw: 75 <math>\mu\text{L}</math></li><li>Total length: 75 mm (3")</li><li>Stem diameter: 2.2 mm (0.1")</li></ul>

#### 100 $\mu\text{L}$ (0.10 mL) capacity




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
139114	Non-sterile	500 pcs/bag, 10 bags/case	5000	<ul style="list-style-type: none"><li>Capacity: 100 <math>\mu\text{L}</math></li><li>Bulb draw: 100 <math>\mu\text{L}</math></li><li>Total length: 83 mm (3.3")</li><li>Stem diameter: 2.3 mm (0.1")</li></ul>

#### 200 $\mu\text{L}$ (0.20 mL) capacity



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
139116	Non-sterile	500 pcs/bag, 10 bags/case	5000	<ul style="list-style-type: none"><li>Capacity: 200 <math>\mu\text{L}</math></li><li>Bulb draw: 200 <math>\mu\text{L}</math></li><li>Total length: 115 mm (4.5")</li><li>Stem diameter: 2.5 mm (0.11")</li></ul>

#### 250 $\mu\text{L}$ (0.25 mL) capacity



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
139118	Non-sterile	500 pcs/bag, 10 bags/case	5000	<ul style="list-style-type: none"><li>Capacity: 250 <math>\mu\text{L}</math></li><li>Bulb draw: 250 <math>\mu\text{L}</math></li><li>Total length: 104 mm (4.1")</li><li>Stem diameter: 2.8 mm (0.11")</li></ul>