

PathAguard ABLiS® Basin Carrier

INSTRUCTIONS FOR USE



1
Before placing ABLiS in ABLiS Basin Carrier make certain that the wheel locking mechanism is applied



2
Make sure that correct height adjustment is completed before moving



3
Holding ABLiS firmly at both ends, place ABLiS in ABLiS Basin Carrier.



4
Be cautious when wheeling carrier, avoid collision with raised objects (e.g. cables/wires)



5
Fill ABLiS with liquid to required level.

PRODUCT CODE

CMCBC10001A: PathAguard ABLiS Basin Carrier



CE