


Mobile Holder for Image Receptors

Cat# MCH-39C



Tall Format Lateral (complete mobile unit) Holds 14" x 36" or 51" image receptors. Fixed position only (non-articulating). Mobile unit includes a vertical range from bottom of image receptor: floor to 54". Posi-lock height adjustable arm with spring loaded pin.

Features 60" chromed steel post. 24" H-Base with locking casters is broad to prevent tipping. Platinum baked enamel finish. Easy to clean base is perfect for OR settings. Improved bilateral clamp safely secures imaging receptors in place.



- Manufactured in Nashville, Tennessee USA
- 14 Day Lead Time
- One Year Warranty
- * Ships in 3 cartons: Pole 12lbs, base 23 lbs , & arm 19 lbs.

