

35X PROFlexx®



70 Weil Way
Wilmington, OH 45177 USA

P 877.733.0911 / 937.382.1451
F 888.388.1349 / 937.382.1191

www.FernoEMS.com

STANDARD FEATURES

- U-Tread Wheels
- 6" Wheels
- Reflective Labels
- Easy-Pull Release Handle
- Integrated Box Beam
- Strengthened C- Channel
- One Handed Side Arm Release
- 1.75" Tube Diameter
- Ratchet Bar Guard
- Dual Wheel Lock
- Gas Assist Backrest
- Powder Coated Frame
- 2 Color Choices: Rescue Red or Electric Blue
- Anodized Patient Deck
- Tested to Industry Standards
- All-Metal Frame
- Built-in Shock Frame
- Stat Trac or 175 Fastener Compatible
- LBS Ready
- One Handed Drop Frame
- 5 Position Drop Frame
- One Handed Safety Bail Release
- Restraints and Mattress Included
- Telescoping Side Lift Handles
- Backrest Adjusts 0 – 65 degrees
- Wider Bed Surface and Head Frame
- Locked-on Mattress
- Foot End Lead Handle

SPECIFICATIONS

Length Max.	79 in. / 2010 mm.
Length Min.	64 in. / 1620 mm.
Width	24 in. / 610 mm.
Height Max.	40 in. / 1020 mm.
Height Min.	14 in. / 350 mm.
Loading Height 1	27 in. / 690 mm.
Loading Height 2	31 in. / 790 mm.
Loading Height 3	35 in. / 890 mm.
Weight	92 lbs. / 42 kg.
Load Limit Loading	700 lbs. / 317 kg.
Load Limit Lowest	1100 lbs. / 500 kg
Shipping Weight	135 lbs. / 61 kg

Ferno reserves the right to change specifications without notice.

ACCESSORIES

CA6099	LBS JR. Bariatric Board
#0314050	LBS JR. Mattress
CA5250	Equipment Bracket
CA5130	Model 513 IV Pole
#0822121	Bariatric Tow Cables
CA2114	Dialysis Kickstand
CA2740	PacRac™
#0314028	Proflexx® Mattress with Transfer Board Pocket
#0313915	417-1 Restraints (Black)
#0313892	430 Restraints 5" (Black)
#0313911	430 Restraints 7" (Black)

OPTIONAL FEATURES

- Universal Sidearms



