

# Model 35-A Cot

 **FERNO**<sup>®</sup>  
When It's Critical<sup>®</sup>



---

70 Weil Way  
Wilmington, OH 45177 USA

P 877.733.0911 / 937.382.1451  
F 888.388.1349 / 937.382.1191

[www.FernoEMS.com](http://www.FernoEMS.com)

---

## STANDARD FEATURES

---

- X-Frame Loading Style
- Fold Down Drop Frame
- Foot End Lead Handle
- Easy-Pull Release Handle
- Side Release Handle
- Lubrication Free
- Easy Roll Wheels
- 8 Bed Positions
- Two Loading Heights
- Backrest Adjusts 0-75 Degrees
- Gas Assist Backrest
- One Hand Swing Down Sidearms
- Safety U bar and Hook
- Tubular Aluminum Construction
- Tested to Industry Standards
- Lightweight
- Built-in Shock Frame
- 175 Fastener Compatible

## SPECIFICATIONS

---

Length Max.	79 in. / 2000 mm.
Length Min.	63 in. / 1600 mm.
Width	24 in. / 610 mm.
Height Max.	40 in. / 1010 mm.
Height Min.	9 in. / 230 mm.
Loading Height 1	28 in. / 710 mm.
Loading Height 2	32 in. / 800 mm.
Weight	74 lbs. / 34 kg.
Load Limit Loading	500 lbs. / 227 kg.
Shipping Weight	105 lbs. / 43 kg.

---

Ferno reserves the right to change specifications without notice.

## ACCESSORIES

---

#0313692	360 Flat Mattress (Burgundy)
#0313691	360 Flat Mattress (Orange)
#0313063	460 Mattress Bolstered (Burgundy)
#0312720	460 Mattress Bolstered (Orange)
#0087156	Model 513 IV Pole
#0818933	PacRac™
#0313915	417-1 Restraints (Black)
#0313914	417-1 Restraints (Orange)
#0313896	417-1 Restraints (Burgundy)
#0313892	430 Restraints 5" (Black)
#0313891	430 Restraints 5" (Orange)
#0313890	430 Restraints 5" (Burgundy)

---



